Modelo BBM-200

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Manual de Instrucciones

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Model BBM-200

Manual Instruction

E_Health Body Slimming & Massage Belt Instruction

Thanks for your purchase of BBM-200 e_health body slimming & massage belt .Before use, please read the specification thoroughly, especially the caution part .so that you will operate it correctly, which should be well keep for your reference at any time .

Contents

Parts Identification

Features of product

Functions

Usage

Cautions

Maintenance and Storage

Specification

Parts Identification

- | | | | |
|----------------------|-----------------------------------|------------------------------|---------------------|
| 1. Massage belt | A. Controller connection wires | | |
| 2. Controller | B. Massage belt connection wires; | C. Adaptor connector; | |
| | D. Power; | E. Auto/Manual; | F. Forward/Reverse; |
| | G. High; | H. Low; | I. Led |
| 3. Outstretched belt | | | |
| 4. DC Adaptor | J. Power plug; | K. Controller connection pin | |

Features of product

- Appearance novelty, use widely.
- Man-style design, simple operation, more well-suited.
- Small volume, different from the former fat reduce products.
- Multi-functions, strong revolving vibration can speed up local fat sport, burning local fat.
- Traditional medicine acupuncture-point therapy combines with modern electronic technology realize local high speed to burnt fat function.
- Three work modes, five level speed adjust, which can annihilate the pertinacity fat perfectly.
- Set Forward/Reverse function can change former product's limitation, bestow local fat repeat to reduce fat properly.
- Eradicate fat and massage to the arm, thigh, calf directly, and the outstretched belt function to lower abdomen, abdomen, waist.
- Body slimming massage belt high-speed circumvolve vibration, make the perfect combination for fat-reduce and massager.
- Time preset function, convenience use.

Functions

With the better living .unreasonable food structure and surplus nutrition .modern people become fat days after days. Recording, there are two-thirds fat people who have to reduce the fat. It is the truth...

The fat that caused by food structure and short for exercises .e_health body slimming massage belt not only in effective allusion to fat .mass fattiness what gluttonize. Crapulence selection partial ion and short for the exercises people grow up bad habit .but build up fit body and health protection massager.

E_health body slimming massage belt high-speed circumvolve vibration can speed up local sport .the wave can stimulate the surplus fat and make surplus fat stored in the body decompose enzymes, speed up fat eradication. Produce heat by stimulate cells movement and immediately dissolve the orange pellicle fatty membrane .meanwhile stimulate blood circulation and lymphatic toxicant elimination ,transmit the dissolved fat to the lymph tube at the same time and discharge out of the body through the excretory system .it can efficiently improve the flexibility of fiber issue ,smooth the skin membrane secretion .increase the body capability of consuming calorie ,curb the growth of new surplus fat and consequently achieve the goal of slimming body and shaping a perfect figure .so using it can increase defecation metes of drainage system meanwhile feeling heat and scratchy .Insist on using it half an hour ,the bestow pot will feel a little aching .there will have a distinct effect .

E_health body slimming massage belt can get effect instantly to aching of waist, leg arthritis and tired body .To have a strong healthy figure using for a long period.

Usage

Please check whether the voltage corresponds with the specification indicated or not before use.

1. Tie up the slimming massage belt on the proper position for body, such as waist, abdomen etc. Tension and adjust to comfortable. If the length for the slimming massage belt is not enough, you can use the outstretched belt.
2. Connect the controller to body slimming massage belt first, than connect to the adaptor, finally

- insert the adaptor plug to the power socket.
3. Press Power switch, enter into work mode 1, the work intensity is 1 level and LED display red. Press HIGH button to increase work intensity, which will be grade increase for each press, the key-press buzzer will stop work when increase to 5 levels. Press LOW button to reduce work intensity, which will be descended for each press, the key-press buzzer will stop work when reduce to 1 level. Hear the FORWARD/REVERSE can not be used.
 4. Under work, when press AUTO/MANUAL button, the work mode will circle among 1, 2 and 3 for press one time.
 5. When the LED display red at work mode 1, FROWARD/REVERSE can not operate. When LED display yellow at work mode 2, FORWARD/REVERSE can not operate too. But when LED display green at work 3, FORWARD/REVERSE can operate.
 6. Press FROWARD/REVERSE button can switch forward/reverse vibration when work mode 3, LED display green.
 7. Massage belt preset the work time is 10 minutes, it will stop operation automatically when the preset time is up.
 8. It is advisable to drink 500cc of lukewarm water before or after use, so as to eliminate the toxic wastes.
 9. When the apparatus is not used, please pull out the plug from the electrical socket.

Cautions

To keep this apparatus the best status, please notice the following problems. Use this massager only for its intended use as described in this manual.

1. Low intensity is suggested at the beginning of use, then improve the intensity gradually, to avoid the discomfort.
2. Please turn off power and pull out the plug after use and before clean it.
3. It only use for fat-reduce, massage and health- protection, please not for others.
4. Use it for 30 minutes at least everyday in order to get fatter-reduce effective.
5. The bestow position will feel heat and scratchy when use it, and there will be a little aching on the second day. They are the natural reaction.
6. Please adjust the speed according to your body condition and needs.
7. The bad habit of Gluttonize, crapulence, selection & partial food are the natural enemy for fat, so when use this apparatus, please notice your bite and sup habit to the perfect effects.
8. To get more effect, please insist on using it everyday.
9. Please stop it if you feel any discomfort or appear the abnormality of apparatus.
10. Do not use it on injured skin.
11. If you are in a medical treatment or involved in any of the following conditions, please follow suggestion of instruction of your doctor:
 - Patient of tumor
 - Patient of chronic disease or serious disease
 - Serious heart disease
 - Being equipped with cardiac pacemaker
 - Being equipped with other metal apparatus as artificial limb
 - the insane
 - the pregnant women
12. Please do not sit or stand on the body slimming massage belt heavily.
13. Do not cast or throw this apparatus.
14. Do not drag its connection wire and power line heavily.
15. Please place this apparatus in a place out of children's reach.
16. With disabled or not clear-headed status person, an assistant shall be in place before use.
17. An apparatus should never be left unattended when plug in, especially if children are present .please take plugs out of the power socket without use.

18. Do not use it in water or damp areas, turn off the power and take its plugs out of socket as soon as drop into the water, do not take out the body slimming massage belt out of the water directly.
19. Do not use near to an oven or any other heating element.
20. Please do not use the fittings and accessory which are not provided by our company.
21. If you take apart the apparatus or replace the fitting in private, we will not guarantee to keep it in good repair, and stop the relative after-sale services.
22. Do not open its slide fastener to keep your security and its longevity. All pose result as to you not obey .our company will be not pay for any duty.
23. If it appears any failure and abnormality, please connect with agent or company which you purchased.
24. Please do not use during the first hour before and after meals.
25. Please drink 500cc of warm water after use.

Caution: the power of this body slimming massage belt is DC24V which reduce by DC adaptor, so it is absolute affective and safe.

Maintenance and storage

1. Do not use benzene, thinner or any other solvent on the massager or fittings.
2. Use a dry cloth when cleaning the massage belt or around the power switch. Never use a damp cloth to clean the massage belt.
3. The body slimming massage belt is specially designed and constructed to be maintenance free. No periodic lubrication or servicing is required.
4. When not in use for a long period of time, it is advisable that the apparatus and its fittings stored in a fire and moisture free environment.
5. Do not store the body slimming massager belt near heat or open flame. Do not leave the massager belt in direct sunlight for a long period of time.

Note: All servicing of this massage belt must be performed by authorized company service person only.

Specification

Description: e_health body slimming & massage belt

Type: BBM-200

Input Voltage: AC100-220V 50/60HZ

Output Voltage: DC24V 2A

Power Consumption: 40W